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Dental Radiographic Examinations: Recommendations for Patient Selection and Limiting Radiation Exposure  
歯科放射線検査：患者選択と放射線被曝制限の推奨事項

## HAND-HELD X-RAY UNITS

Hand-held, battery-powered x-ray systems are available for intra-oral radiographic imaging. The hand-held exposure device is activated by a trigger on the handle of the device. However, dosimetry studies indicate that these hand-held devices present no greater radiation risk than standard dental radiographic units to the patient or the operator. No additional radiation protection precautions are needed when the device is used according to the manufacturer's instructions. These include: 1. holding the device at mid-torso height, 2. orienting the shielding ring properly with respect to the operator, and 3. keeping the cone as close to the patient's face as practical. If the hand-held device is operated without the ring shield in place, it is recommended that the operator wear a lead apron.

All operators of hand-held units should be instructed on their proper storage. Due to the portable nature of these devices, they should be secured properly when not in use to prevent accidental damage, theft, or operation by an unauthorized user. Hand-held units should be stored in locked cabinets, locked storage rooms, or locked work areas when not under the direct supervision of an individual authorized to use them. Units with user-removable batteries should be stored with the batteries removed. Records listing the names of approved individuals who are granted access and use privileges should be prepared and kept current.

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